DSM Food Specialties B.V.

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Delvotest® BLF

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DSM Food Specialties B.V.

Delvotest® BLF

1. Take the strip container out of refrigeration (4-8°C recommended storage conditions) and allow the container to warm up to room temperature

2. Detach one or more ampoules. Take care that the aluminum foil from the remaining tests is not damaged. Perforate the ampoules carefully.

3. Take a fresh disposable pipette for each milk sample. **Do not re-use pipettes.** Do not touch the tip-end, which will be in contact with the milk.

4. Add the milk into the pipette by squeezing the smaller upper bulb **once**, hold it, dip the pipette tip about 1 cm into the milk sample. Then release pressure on the bulb and the pipette (stem) will fill itself with the appropriate volume (0.15 ml) of milk.

Note: After pipetting, the little surplus of milk is caught into the pipette reservoir (the small lower bulb). Repeat the pipetting when there is no milk in the reservoir.

5. Transfer the milk samples by gently and totally squeezing the same upper bulb, adding the milk straight onto the agar medium. **Note:** The little surplus of milk in the lower bulb will remain.

Shake for 10 seconds

6. Check the temperature of the incubator (64°C +/- 2°C). Put the ampoules into the incubator.

Wait 2 minutes

Shake for 5 seconds


7. Take one test strip; immediately close the container containing the strips. Hold the test strip from the upper side (green area). Insert one test strip per ampoule.

8. Remove the strip after **3 minutes** and read the result.

Control < Test → Negative

Control ≥ Test → Positive
(It contains antibiotics)

Points of Attention:

- 1- As this test is extremely sensitive to antibacterial substances (β-lactams) any contamination with these substances should be prevented. It is advisable to thoroughly wash and dry hands before starting the test procedure. Use a clean table.
- 2- Do not handle the test roughly since this may damage the test strips. This may also affect the quality of the test during reading of the test results.
- 3- After removing the test strips from their container, immediately close the container and return it to the recommended storage conditions to avoid possible effects from humidity which shortens the strip shelf life.

Although diligent care has been used to ensure that the information provided herein is accurate, nothing contained herein can be construed to imply any representation or warranty for which we assume legal responsibility, including without limitation any warranties as to the accuracy, currency or completeness of this information or of non-infringement of third party intellectual property rights. The content of this document is subject to change without further notice. Please contact us for the latest version of this document or for further information. Since the user's product formulations, specific use applications and conditions of use are beyond our control, we make no warranty or representation regarding the results which may be obtained by the user. It shall be the responsibility of the user to determine the suitability of our products for the user's specific purposes and the legal status for the user's intended use of our products.

Sensitivity list


DSM Food Specialties B.V.

Delvotest® BLF

The Delvotest® BLF is a fast test for the detection of β -lactam antibiotics. The Delvotest® BLF is therefore sensitive to all drug containing a β -lactam ring (see figure), including

- penams (pencillins) extended and narrow spectrum,
- penams β -lactamase resistant and β -lactamase sensitive;
- cephalosporines of 1st, 2nd, 3rd, 4th, and 5th generation.

The list below is not exhaustive. The sensitivity of the Delvotest BLF is given only for the most commonly used antibiotics.

Basic structure of β -lactam drugs detected by the Delvotest® BLF

Family	Molecules	MRL EU	Delvotest BLF
Penams (Penicillins)	amoxicillin	4	3
	ampicillin	4	5
	cloxacillin	30	17
	dicloxacillin	30	24
	nafcillin	30	54
	oxacillin	30	36
	penethamate	4	*
	penicillin G	4	3
Cephalosporines (1st, 2nd, 3rd, 4th, 5th generation)	cefacertrile	125	*
	cefalexin	100	34
	cefalonium	20	3
	cefapirin	60	*
	cefazolin	50	10
	cefoperazone	50	5
	cefquinome	20	11
	ceftiofur [#]	100	4
	cefuroxime	?	*

*detected but exact sensitivity not known

[#]ceftiofur pur (parent drug)

Although diligent care has been used to ensure that the information provided herein is accurate, nothing contained herein can be construed to imply any representation or warranty for which we assume legal responsibility, including without limitation any warranties as to the accuracy, currency or completeness of this information or of non-infringement of third party intellectual property rights. The content of this document is subject to change without further notice. Please contact us for the latest version of this document or for further information. Since the user's product formulations, specific use applications and conditions of use are beyond our control, we make no warranty or representation regarding the results which may be obtained by the user. It shall be the responsibility of the user to determine the suitability of our products for the user's specific purposes and the legal status for the user's intended use of our products.